


Aspicath Thrombus Aspiration Catheter



*Biometrix Thrombus Aspiration Catheters
provide the best solution for thrombus removal*

Aspicath Thrombus Aspiration Catheter

The Biometrix Thrombus Aspiration Catheter incorporates a smooth tip design to increase pushability and a unique hydrophilic coating to escalate trackability. The Aspicath's excellent pushability is also due to the innovative shaft design and the kink resistance, specifically designed for aspiration catheters.



Aspiration Catheter Set

Code	Description
BG-A6F-ASPICATH	6F Aspiration Catheter
BG-A7F-ASPICATH	7F Aspiration Catheter

Qty/Box: 1

Technical Data	Aspicath 6F	Aspicath 7F
Range	6F	7F
Wire compatibility	0.014" (0.36mm)	0.014" (0.36mm)
Working length	145cm	145cm
Distal part flexibility/trackability	Pebax flexibility	Pebax flexibility
Proximal part push/kinking resistance	Reinforced PEEK for support	Reinforced PEEK for support
Outer diameter (O.D.)	0.067" (1.70mm)	0.078" (1.98mm)
Extraction lumen area (Proximal)	0.95mm ²	1.54mm ²
Extraction lumen area (Distal)	0.93mm ²	1.43mm ²
Aspiration speed	1.6cc/s.	2.8cc/s.

