

ACTISEAL



bipolar fingerswitch™




- Fingerswitch controlled – no footswitch
- Non-stick sealing tips
- Perfectly balanced for fast decisive surgery



control without **compromise**



bipolar fingerswitch™




Ordering details:

- Actiseal with 0.8mm tips, 4 metre cable and fingerswitch control - **PMS8208**
- Actiseal with 1.2mm tips, 4 metre cable and fingerswitch control - **PMS8212**

Sterile packed in 20's

- Tips enable tissue sealing by conducting the heat responsible for tissue adhesion away from the operative site
- Virtually no tissue adhesion

For use with Valleylab, Megadyne, Conmed and Bovie electro-surgery generators



3 pin plug for Valleylab, Megadyne and Conmed electro-surgery generators

Actiseal with 1.2mm tip illustrated



International patents, design registrations and trade marks applied for on both product and tubing

Tel: +44 (0)1453 752626 Email: sales@prima-medical.com Web: www.prima-medical.com

Prima Medical Limited, Mill 16, Bath Road Trading Estate, Stroud, Gloucestershire, GL5 3QF, United Kingdom



Prima
Inspiring innovation