


adjust

adjustable bipolar forceps

- 34 increments in 1mm stops
- Non-stick coating, minimises eschar build up
- Unique design allows the surgeon to have multiple forceps lengths
- Pressure switching, allowing accurate safe seal
- Fingerswitch activation
- Single use





adjust

adjustable bipolar forceps

- Heat conducting nonstick tips
- Instantly has a range of forceps to hand reducing the number of forceps required during surgery
- Product options to suit you:
 - Option 1: with switch on forceps
 - Option 2: without switch; to be used with footswitch

Finger grips move in 34 increments in 1mm stops



Ordering details:

- 1.0mm tips - PMS1500 - boxed 10
- 1.5mm tips - PMS1600 - boxed 10
- 2.0mm tips - PMS1700 - boxed 10
- Adjust nonstick forceps
- 3m cable with standard 3 pin connector