


Ceto

RX PTCA Dilatation Balloon Catheter



Biometrix RX PTCA dilatation balloon catheter provides the best solution for coronary angioplasty procedures

The image features a blue gradient background with a white horizontal band at the top. A thin, dark, flexible catheter is shown in a large, elegant loop that spans across the middle and right side of the frame. The catheter's tip is visible on the right, showing a tapered, needle-like end. The overall aesthetic is clean and professional, typical of medical product marketing.

Ceto

RX PTCA Dilatation Balloon Catheter


Biometrix RX PTCA Dilatation balloon catheter is semi-compliant. It is used for dilatation of clogged coronary blood vessels in angioplasty procedures. The catheter has excellent pushability features and a flexible tapered tip. With optimal trackability and an innovative hydrophilic coating (FC) it is fit for the tightest occlusion crossings.


		Diameter in mm											
		1.50	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00	4.50	5.00
Length in mm	8 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	10 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	15 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	20 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	25 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	30 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	40 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

✓ Available

		Diameter in mm											
		1,50	2,00	2,25	2,50	2,75	3,00	3,25	3,50	3,75	4,00	4,50	5,00
Pressure in ATM	2	1,39	1,85	2,09	2,32	2,55	2,78	3,01	3,24	3,48	3,71	4,17	4,63
	4	1,44	1,92	2,16	2,40	2,64	2,88	3,12	3,36	3,60	3,84	4,32	4,80
	6	1,50	2,00	2,25	2,50	2,75	3,00	3,25	3,50	3,75	4,00	4,50	5,00
	8	1,54	2,05	2,31	2,56	2,82	3,07	3,33	3,58	3,84	4,09	4,60	5,12
	10	1,57	2,09	2,35	2,61	2,87	3,13	3,39	3,65	3,91	4,17	4,69	5,22
	12	1,59	2,12	2,39	2,65	2,92	3,18	3,45	3,71	3,98	4,24	4,77	5,30
	14	1,62	2,16	2,43	2,70	2,97	3,24	3,51	3,78	4,05	4,32	4,86	5,40
	16	1,65	2,20	2,48	2,75	3,03	3,30	3,58	3,85	4,13	4,40	4,95	5,50
	18	1,68	2,24	2,52	2,80	3,08	3,36	3,64	3,92	4,20	4,48	5,04	5,60

Nominal Pressure (NP)
Rated Burst Pressure (RBP)
Balloon size measured at 37°C





Innovative flexible tip for the tightest occlusion crossing

ORDERING INFORMATION

Code	Balloon Diameter	Balloon Length
BG-C15010-CETO	1,5mm	10mm
BG-C15015-CETO	1,5mm	15mm
BG-C15020-CETO	1,5mm	20mm
BG-C15030-CETO	1,5mm	30mm
BG-C20010-CETO	2mm	10mm
BG-C20015-CETO	2mm	15mm
BG-C20020-CETO	2mm	20mm
BG-C20030-CETO	2mm	30mm
BG-C22510-CETO	2,25mm	10mm
BG-C22515-CETO	2,25mm	15mm
BG-C22520-CETO	2,25mm	20mm
BG-C22530-CETO	2,25mm	30mm
BG-C25010-CETO	2,5mm	10mm
BG-C25015-CETO	2,5mm	15mm
BG-C25020-CETO	2,5mm	20mm
BG-C25030-CETO	2,5mm	30mm
BG-C30010-CETO	3mm	10mm
BG-C30015-CETO	3mm	15mm
BG-C30020-CETO	3mm	20mm
BG-C30030-CETO	3mm	30mm
BG-C35010-CETO	3,5mm	10mm
BG-C35015-CETO	3,5mm	15mm
BG-C35020-CETO	3,5mm	20mm
BG-C35030-CETO	3,5mm	30mm
BG-C40010-CETO	4mm	10mm
BG-C40015-CETO	4mm	15mm
BG-C40020-CETO	4mm	20mm
BG-C40030-CETO	4mm	30mm

Qty/Box: 1