


DRAINAGE & NEPHROSTOMY CATHETERS

DRAINAGE CATHETER SETS Hydrophilic

Pigtail Catheter Sets

One Step Pigtail with Safety String Lock



Code	French Size	Length (In Open Position)	Accepts Guide Wire Diameter	Flexible Stiffening Cannula
NC-06191	6F	20cm	0.035"	-
NC-06241	6F	25cm	0.035"	-
NC-06291	6F	30cm	0.035"	-
NC-06391	6F	40cm	0.035"	-
NC-08191	8F	20cm	0.038"	-
NC-08241	8F	25cm	0.038"	-
NC-08291	8F	30cm	0.038"	-
NC-08391	8F	40cm	0.038"	-
NC-10191	10F	20cm	0.038"	-
NC-10241	10F	25cm	0.038"	-
NC-10291	10F	30cm	0.038"	-
NC-10391	10F	40cm	0.038"	-
NC-12191	12F	20cm	0.038"	-
NC-12241	12F	25cm	0.038"	-
NC-12291	12F	30cm	0.038"	-
NC-12391	12F	40cm	0.038"	-
NC-14191	14F	20cm	0.038"	-
NC-14241	14F	25cm	0.038"	-
NC-14291	14F	30cm	0.038"	-
NC-14391	14F	40cm	0.038"	-
NC-08192	8F	20cm	0.038"	+
NC-08242	8F	25cm	0.038"	+
NC-08292	8F	30cm	0.038"	+
NC-08392	8F	40cm	0.038"	+
NC-10192	10F	20cm	0.038"	+
NC-10242	10F	25cm	0.038"	+
NC-10292	10F	30cm	0.038"	+
NC-10392	10F	40cm	0.038"	+
NC-12192	12F	20cm	0.038"	+
NC-12242	12F	25cm	0.038"	+
NC-12292	12F	30cm	0.038"	+
NC-12392	12F	40cm	0.038"	+
NC-14192	14F	20cm	0.038"	+
NC-14242	14F	25cm	0.038"	+
NC-14292	14F	30cm	0.038"	+
NC-14392	14F	40cm	0.038"	+

Set Contents:

- Drainage Catheter
- Safety String Lock
- Metal Stiffening Cannula
- Trocar Needle
- Split Straightener
- *Flexible Stiffening Cannula
(only for 8-14Fr catheter)

Customized sets available with the following optional contents:

- Connecting tube with 3 way stopcock
- Scalpel
- Drainage bag
- Securement


Qty/Box: 5

DRAINAGE & NEPHROSTOMY CATHETERS

DRAINAGE CATHETER SETS Hydrophilic

Pigtail Catheter Sets

Closed Loop Pigtail with Safety String Lock (Seldinger Type)



Code	French Size	Length (In Open Position)	Accepts Guide Wire Diameter	Flexible Stiffening Cannula
NC-06194	6F	20cm	0.035"	-
NC-06244	6F	25cm	0.035"	-
NC-06294	6F	30cm	0.035"	-
NC-06344	6F	35cm	0.035"	-
NC-06394	6F	40cm	0.035"	-
NC-08194	8F	20cm	0.038"	+
NC-08244	8F	25cm	0.038"	+
NC-08294	8F	30cm	0.038"	+
NC-08344	8F	35cm	0.038"	+
NC-08394	8F	40cm	0.038"	+
NC-10194	10F	20cm	0.038"	+
NC-10244	10F	25cm	0.038"	+
NC-10294	10F	30cm	0.038"	+
NC-10344	10F	35cm	0.038"	+
NC-10394	10F	40cm	0.038"	+
NC-12194	12F	20cm	0.038"	+
NC-12244	12F	25cm	0.038"	+
NC-12294	12F	30cm	0.038"	+
NC-12344	12F	35cm	0.038"	+
NC-12394	12F	40cm	0.038"	+
NC-14194	14F	20cm	0.038"	+
NC-14244	14F	25cm	0.038"	+
NC-14294	14F	30cm	0.038"	+
NC-14344	14F	35cm	0.038"	+
NC-14394	14F	40cm	0.038"	+

Set Contents:

- Closed Loop Pigtail Drainage Catheter
- Safety String Lock
- Metal Stiffening Cannula
- Split Straightener
- *Flexible Stiffening Cannula (only for 8-14Fr catheter)

Customized sets available with the following optional contents:


- Connecting tube with 3 way stopcock
- Scalpel
- Drainage bag
- Securement

Qty/Box: 5

DRAINAGE CATHETER SETS Hydrophilic

Pigtail Catheter Sets

Closed Loop Pigtail (One Step Type)



Code	French Size	Length (In Open Position)	Accepts Guide Wire Diameter	Safety String Lock
NC-06195	6F	20cm	0.035"	+
NC-06245	6F	25cm	0.035"	+
NC-06295	6F	30cm	0.035"	+
NC-06395	6F	40cm	0.035"	+
NC-08195	8F	20cm	0.038"	+
NC-08245	8F	25cm	0.038"	+
NC-08295	8F	30cm	0.038"	+
NC-08395	8F	40cm	0.038"	+
NC-10195	10F	20cm	0.038"	+
NC-10245	10F	25cm	0.038"	+
NC-10295	10F	30cm	0.038"	+
NC-10395	10F	40cm	0.038"	+
NC-12195	12F	20cm	0.038"	+
NC-12245	12F	25cm	0.038"	+
NC-12295	12F	30cm	0.038"	+
NC-12395	12F	40cm	0.038"	+
NC-14195	14F	20cm	0.038"	+
NC-14245	14F	25cm	0.038"	+
NC-14295	14F	30cm	0.038"	+
NC-14395	14F	40cm	0.038"	+
NC-06196	6F	20cm	0.035"	-
NC-06246	6F	25cm	0.035"	-
NC-06296	6F	30cm	0.035"	-
NC-06396	6F	40cm	0.035"	-
NC-08196	8F	20cm	0.038"	-
NC-08246	8F	25cm	0.038"	-
NC-08296	8F	30cm	0.038"	-
NC-08396	8F	40cm	0.038"	-
NC-10196	10F	20cm	0.038"	-
NC-10246	10F	25cm	0.038"	-
NC-10296	10F	30cm	0.038"	-
NC-10396	10F	40cm	0.038"	-
NC-12196	12F	20cm	0.038"	-
NC-12246	12F	25cm	0.038"	-
NC-12296	12F	30cm	0.038"	-
NC-12396	12F	40cm	0.038"	-
NC-14196	14F	20cm	0.038"	-
NC-14246	14F	25cm	0.038"	-
NC-14296	14F	30cm	0.038"	-
NC-14396	14F	40cm	0.038"	-

Set Contents:

- Closed Loop Pigtail Drainage Catheter
- Metal Stiffening Cannula
- Trocar Needle
- Split Straightener

Customized sets available with the following optional contents:

- Connecting tube with 3 way stopcock
- Scalpel
- Drainage bag
- Securement


Qty/Box: 5

DRAINAGE & NEPHROSTOMY CATHETERS

DRAINAGE CATHETER SETS Hydrophilic

Biliary Catheter Sets

Biliary Catheter with Radiopaque Band and Safety String Lock



Code	French Size	Length (In Open Position)	Accepts Guide Wire Diameter	Flexible Stiffening Cannula
NC-0619B	6F	20cm	0.035"	-
NC-0624B	6F	25cm	0.035"	-
NC-0629B	6F	30cm	0.035"	-
NC-0639B	6F	40cm	0.035"	-
NC-0644B	6F	45cm	0.035"	-
NC-0819B	8F	20cm	0.038"	+
NC-0824B	8F	25cm	0.038"	+
NC-0829B	8F	30cm	0.038"	+
NC-0839B	8F	40cm	0.038"	+
NC-0844B	8F	45cm	0.038"	+
NC-1019B	10F	20cm	0.038"	+
NC-1024B	10F	25cm	0.038"	+
NC-1029B	10F	30cm	0.038"	+
NC-1039B	10F	40cm	0.038"	+
NC-1044B	10F	45cm	0.038"	+
NC-1219B	12F	20cm	0.038"	+
NC-1224B	12F	25cm	0.038"	+
NC-1229B	12F	30cm	0.038"	+
NC-1239B	12F	40cm	0.038"	+
NC-1244B	12F	45cm	0.038"	+
NC-1419B	14F	20cm	0.038"	+
NC-1424B	14F	25cm	0.038"	+
NC-1429B	14F	30cm	0.038"	+
NC-1439B	14F	40cm	0.038"	+
NC-1444B	14F	45cm	0.038"	+

Set Contents:

- Biliary Drainage Catheter with Radiopaque Band
- Safety String Lock
- Metal Stiffening Cannula
- Split Straightener
- MLL Injection Site
- *Flexible Stiffening Cannula (only for 8-14Fr catheter)


Customized sets available with the following optional contents:

- Connecting tube with 3 way stopcock
- Scalpel
- Drainage bag
- Securement

Qty/Box: 5

DRAINAGE CATHETER SETS


Drainage Catheter Insertion Set



Code	Description	Qty/Box
HD-N6410	Drainage Catheter Insertion Set	5

Set Contents:

- 21Gax15cm Access Needle • OD: 0.5 mm Stylet • 0.018"x 65cm Special Kink Resistant Guidewire
- 22 Gauge Stiffening Cannula • 4Fr x 20cm Dilator • 6Fx18cm Sheath with Hydrophilic Coating



Code	Description	Qty/Box
HD-N6411	Drainage Catheter Insertion Set with 0.035"x150cm Amplatz Guidewire	5
HD-N6412	Drainage Catheter Insertion Set with 0.035"x80cm Amplatz Guidewire	5

Set Contents:

- 21Gax15cm Access Needle • OD: 0.5 mm Stylet • 0.018"x 65cm Special Kink Resistant Guidewire
- 22 Gauge Stiffening Cannula • 4Fr x 20cm Dilator • 6Fx18cm Sheath with Hydrophilic Coating
- 0.035"x 80/150cm Amplatz Guidewire

10mm from tip Radiopaque Band on Sheath

Drainage Catheter Connectors

Code	Length	Stopcock	Qty/Box
DQ-0305	5cm	-	50
DQ-0325	25cm	-	50
DQ-0430	40cm	+	100

