

MACRO SEAL




monopolar tissue sealing forceps

- Fingerswitch controlled – no footswitch
- Non-stick sealing tips
- Perfectly balanced for fast, decisive surgery



control without **compromise**

MACRO SEAL



Ordering details:


- Macro Seal with 1,2mm tips, 4 metre cable and fingerswitch control - **PMS8412**

Sterile packed in 20's

monopolar tissue sealing forceps

- Tips enable tissue sealing by conducting the heat responsible for tissue adhesion away from the operative site
- **nonstick** finish ensures adhesion free coagulation
- Macro Seal is activated by pressing buttons on the side - continued pressure activates the electro-surgery power

For use with all electro-surgery generators



3 pin plug for all popular electro-surgery



International patents, design registrations and trade marks applied for on both product and tubing

Tel: +44 (0)1453 752626 Email: sales@prima-medical.com Web: www.prima-medical.com

Prima Medical Limited, Mill 16, Bath Road Trading Estate, Stroud, Gloucestershire, GL5 3QF, United Kingdom



Prima
Inspiring innovation