


# SUPER VAC



## telescopic smoke evacuation fingerswitch


- Soft touch switching buttons
- Puts vacuum where you need it
- Handle extends from 20mm to 85mm enabling fingertip and deep cavity surgery to be performed using only one product
- length is adjusted in seconds



**control** without **compromise**



telescopic, high volume smoke evacuation fingerswitch



- Integral smoke evacuation with high flexibility tubing
- Comfortable triangular grip
- Easy push/pull control carriage

Product	Quantity per case
<b>PMS6300</b> SuperVac smoke evacuation fingerswitch with non stick blade electrode	20
<b>PMS6301</b> SuperVac (Plastic) smoke evacuation fingerswitch with tungsten needle electrode	20
<b>PMS6302</b> SuperVac smoke evacuation fingerswitch without electrode	20
<b>PMS6304</b> SuperVac smoke evacuation fingerswitch with SOLO dynamic electrode	20
<b>PMS6305</b> SuperVac smoke evacuation fingerswitch with 4mm non stick blade electrode	20

*All products with 4metre cable and terminated with three pin electrosurgery plug*



International patents, design registrations and trade marks applied for on both product and tubing

+44 (0)1453 752626 sales@prima-medical.com www.prima-medical.com

Prima Medical Limited, Mill 16, Bath Road Trading Estate, Stroud, Gloucestershire, GL5 3QF, United Kingdom



**Prima**  
Innovation in Electrosurgery