


URIMETRIX™ Urine Drainage Program

Closed Urine-Meter System


500ml system connected to
2000ml urine bag with emptying valve.

Product Specifications


- 500ml measuring scale.
- Measuring scale resolution:
0-20ml +/- 1ml
20-110ml +/- 5ml
110-140ml +/- 10ml
140-200ml +/- 20ml
200-500ml +/- 25ml
- Integral emptying tap.
- Overflow segment.
- Anti-reflux valve.
- 2000ml exchange bag
with an emptying valve (UA-0003).
- Built-in sampling port
in the container (UA-0001).



Needle Free Sampling Port



Built-In Sampling Port with Cap
(UA-0001)



2000cc URX Exchange Bag
(UA-0003)

Code	Description
UA-0000	500ml System with Needleless Sampling Port and 150cm Connecting Tube
UA-0001	500ml System with Needleless Sampling Port, Extra Top Sampling Port and 150cm Connecting Tube
UA-0002	2000ml Urimetrix Exchange Bag
UA-0003	2000ml Urimetrix Exchange Bag with Emptying Tap

Qty/Box: 10, UA-0002, UA-0003: 50