

# Intra-Abdominal Pressure (IAP) Monitoring Set

The Intra-Abdominal Pressure Monitoring Set measures intra-abdominal pressure through a disposable transducer by infusing saline into the urine bladder.


A simple and easy-to-use set allowing as many accurate readings as needed from the same patient. The set is connected between the Foley Catheter and the drainage bag or container.

## Intra-abdominal elevated pressure can indicate:

- Increased intrathoracic, intracranial and intracardiac filling pressures
- Decreased left ventricular, chest wall and total respiratory system compliance
- Atelectasis and increased extravascular lung water
- Acute renal failure
- Bacterial translocation and multiple organ system failure

## Set Contents:

- Double Check Valve
- 50ml Syringe
- 3 Way Stopcock
- Large Bore 3 Way Stopcock
- Stepped Adaptor
- Funnel Connector



## Ordering Information

Code	Description
CA-0200	Intra-Abdominal Pressure Monitoring Set
Qty/Box: 20	



Paardeweide 36 P.O.B. 7191, 4800 GD Breda, The Netherlands Tel: 31-76-5422720 Fax: 31-76-5422716  
sales@biometrixmedical.com www.biometrixmedical.com