





# Post Operative Pain Management Simple and effective




CIMPAX<sup>®</sup> Wound Infiltration Catheter Kit  
Complete local anaesthesia kit for post-surgical pain relief  
Including CIMPAX<sup>®</sup> C-CAT catheter


CIMPAX® Wound Infiltration Catheter Kit, is the complete kit for pre-, peri- and post-operative pain management. The kit contains needles for infiltration of the surgical site during the procedure and the C-CAT catheter for distribution of anaesthetics to the surgical site post-operatively. The catheter can be used with infusion pumps as well as for bolus injections. For bolus injections the C-CAT catheter offers minimal resistance due to its unique design.




- Precise, direct and concentrated pain relief
- Reduces or removes the use of narcotics
- Less or no side effects caused by narcotics
- Less pain (improved VAS score)
- Easier patient mobilisation
- Faster recovery, shortens hospital stay
- Increased patient satisfaction




1  
Infiltration of the surgical site



2  
Introduction of the catheter sheath



3  
Introducing catheter through the sheath



4  
Closing the wound and fixation of catheter

## Ordering information

### C-CAT

Catheter with 6 holes covering 2,5 cm Ref. C44-21-005  
 Catheter with 12 holes covering 6,5 cm Ref. C44-21-010  
 Catheter with 30 holes covering 15 cm Ref. C44-21-015  
 Catheter with 42 holes covering 21 cm Ref. C44-21-018  
 Catheter with 60 holes covering 30 cm Ref. C44-21-020  
 Packaging: Box of 14 kits, sterile

### C-CAT, light

Catheter with 6 holes covering 2,5 cm Ref. C44-21-105  
 Catheter with 12 holes covering 6,5 cm Ref. C44-21-110  
 Catheter with 30 holes covering 15 cm Ref. C44-21-115  
 Catheter with 42 holes covering 21 cm Ref. C44-21-118  
 Catheter with 60 holes covering 30 cm Ref. C44-21-120  
 Packaging: Box of 14 kits, sterile