

Single use Insufflation Needle



C-NEEDLE

The CIMPAX[®] C-NEEDLE offers simple and secure pneumoperitoneum for laparoscopic surgery. The single use needle insures cleanness and sterility for every procedure and is a cost effective alternative to reusable needles.

- Luer lock connection with on-off valve
- Spring loaded
- Single use
- Available in 2 sizes, 120 and 150 mm
- Packaging: Box of 10 pcs, sterile

Ref. code 120 mm: C23-12-121

Ref. code 150 mm: C23-12-151

