


Single use Electrosurgical Pencil



C-PEN

The CIMPAX® Single Use Electrosurgical Pencil C-PEN is a high quality pencil. The electrode is securely fitted with a hexagonal connection securing the electrode in the required position. C-PEN can be used with any electrosurgical generator available.

- Comfortable design and light weight
- Colour marked buttons
- Hexagonal electrode connection
- 2.35mm electrode
- 3 pin US-connector
- Single use

C-PEN with 3 meter cable Ref. code C31-08-001

Packaging: Box of 25 pcs, Sterile

C-PEN with 5 meter cable Ref. code C31-08-501

Packaging: Box of 20 pcs, Sterile

CE0543