




Laparoscopic instruments for perfect performance

- Quality
- Precision
- Reliable



Monopolar Hand Instruments
for laparoscopic surgery
Single use

The CIMPAX® Monopolar Laparoscopic Hand Instruments offer everything the professional physician requires to ensure perfect performance through out the complete procedure.



C-CUT 11 or 17 mm blades



C-DISSECT 20 mm jaw length



C-GRASP 18 mm jaw length



C-GRASP 22 mm jaw length



C-CLINCH 22 mm jaw length



C-BAB 22 mm jaw length

- Lightweight well balanced instruments
- Comfortable, easy open and close handle
- Naturally located lock for ratchet feature
- Free and simple rotation of jaws
- Double hinge design ensures easy and precise jaw motion and maximum opening
- For 5 mm ports. Shaft length 330 mm

Ordering information

C-CUT	Laparoscopic Scissors, Metzembaum Curved	17 mm blades	Box of 12, sterile	Ref. C23-11-331
C-CUT	Laparoscopic Scissors, Metzembaum Curved	11 mm blades	Box of 12, sterile	Ref. C23-11-337
C-DISSECT	Laparoscopic Dissector, Maryland Curved	20 mm jaws	Box of 12, sterile	Ref. C23-11-334
C-GRASP	Laparoscopic Grasper, Atraumatic with Ratchet	18 mm jaws	Box of 12, sterile	Ref. C23-11-332
C-GRASP	Laparoscopic Grasper, Atraumatic with Ratchet	22 mm jaws	Box of 12, sterile	Ref. C23-11-336
C-CLINCH	Laparoscopic Grasper, Clinch Grasper with Ratchet	22 mm jaws	Box of 12, sterile	Ref. C23-11-333
C-BAB	Laparoscopic Forceps, Babcock with Ratchet	22 mm jaws	Box of 12, sterile	Ref. C23-11-335