






Electrosurgical pencil with smoke evacuation




CIMPAX® C-VAC TORNADO™
Latex, lead, PVC and phthalate FREE
Single use

CIMPAX® C-VAC TORNADO™ the ultimate combination of electrosurgery and removal of hazardous smoke in one instrument. The pencil has a unique triangular design which combined with the silicone cover offers the surgeon perfect ergonomic and safety grip. The ultra soft micro switch push buttons gives further comfort for the user due the exceptional tactile feedback. Since the C-VAC TORNADO™ is Latex, lead, PVC and phthalate free you can perform surgery without concerns for your health, your colleagues and patient's health or the environment.



Adaptor with 15/22 mm male/female connector and build-inn specimen filter



3-pin connector for HF - generators



Two versions all in one electrode and suction tip extenders with unlimited extension possibilities

- Unique triangular design from tip to end
- Ultra soft push buttons for enhanced comfort
- Maximum safety grip silicone cover
- Transparent tip for clear view of the electrode
- LATEX, LEAD, PVC and PHTHALATE FREE
- Flexible tubing with integrated ESU cable
- Ready to use, no assembling of suction tip or electrode required.

Ordering information

C-VAC TORNADO™ Integrated smoke evacuation pencil. 4 meter lead, PVC and phthalate free ESU cable with US-3-pin connector and 3 meter PVC and phthalate free tubing with combined 22/15/9 mm adaptor. Single use.

With standard blade electrode	Ref. C31-19-001	Extender, short suction tip	Ref. C31-15-145
With coated blade electrode	Ref. C31-19-101	Extender, long suction tip	Ref. C31-15-150
With silicone blade electrode	Ref. C31-19-201	Packaging: Box of 15 pcs., sterile	