




# OBSTETRICS





1652

1651

## Pregdate Calender I 1652/1651

Wing Plast's Pregdate Calender Table Size 1652 and Pocket Size 1651 are pregnancy calendars for indicating the date of arrival of the expected baby. A well-known instrument that every midwife must have to be able to work efficiently.

- Easy to use
- Clear and accurate scale
- IFU on the back of the product



A1390

## Hegenerger Speculum I A1390

Hegenerger Speculum is a sterile, disposable plastic speculum for postpartum sutures. The ergonomic design makes the treatment intuitive and maximises patient comfort. Hegenerger Speculum gives free view of the work area, stays in place by itself which frees hands for suturing and eliminates the need for an assistant to keep the area open. Sterile, in unit package.

- Free view of the suturing area
- Higher quality of suturing
- Low friction surface for easy insertion
- Smooth surfaces and smoothly rounded corners
- Reduces patient discomfort
- Faster suturing time
- No limited posture during work
- Eliminates need for assistant
- Diminished risk of needle stick injuries



## Pregnancy Measuring tape reusable | 1654

Wing Plast's Pregdate Measurer 1653 is used together with the Pregdate calendars for indicating the date of arrival of the expected baby. A well-known instrument that every midwife must have to be able to work efficiently.

- Easy to use
- Clear and accurate scale



## Pregnancy Measuring tape disposable I 1655

Wing Plast's Pregnancy Measuring Tape Disposable is a flexible measuring tape that is easy to handle and store. The measuring tape measures 76 cm. White with black numerals.

- Clear scale
- 20 mm wide
- 760 mm length
- Disposable
- Easy to handle

WM003-01



WM004-01



## Wegaa Symfys-Fundus measuring I WM003-01/WM004-01

Symfys-Fundus measuring tape is a disposable measuring tape for measuring pregnant stomach.

Head measuring tape is a disposable measuring tape for measuring the baby's head circumference.

WEGAA® Symphys-Fundus measuring tape and packaging are made of environmental certified paper. They are packaged in 80 pieces. The disposable measuring tapes are clean and stored in packages that are easy to carry. They are easy to open and convenient to use and clear to read the description. The packaging can be lying or hanging on the wall. WEGAA® disposable measuring tapes are Swedish manufactured and CE marked.

- Disposable
- Swedish-made
- Environmentally certified paper
- User friendly



9610-2

## CTG band Disposable I 9610-2

Wing Plast disposable strap for CTG are 1.3 m long and it is equipped with an extra button. 2-pack with 1 blue and 1 pink ribbon.


- Elastic
- 50 mm wide
- Disposable
- 2 packs



## CTG-band

Wing Plast's CTG band with buttonholes is available in 3 different widths; 32 mm, 50 mm and 100 mm. Multi-purpose tapes that can be washed. Length 15 m.

- Latex Free
- Reusable
- Available in 3 different widths




**30952090**

## **Double button for CTG band 30952090**

Double button for CTG band

- Polystyrene




**421/3**

## **Pozzi pliers I 421/3**

Pozzi pliers take hold of the Cervix





9606

## Gyn table protection I 9606

Collection bag at delivery and gyn-op.  
The bag is graded for estimation by amount of  
amniotic fluids or other such fluids Gynae-Op.

Well-sized size for practical use and handling.  
Made of transparent PE plastic.

- Dimensions, graded
- Very appreciated





9601

## CetroHook I 9601

Cetro Hook made of environmentally friendly plastic that provides a firm grip when handling. Easy to use, with protected tip that does not harm the mother or baby. With the right length for good and safe use.

- Sterile
- Rounded top
- Sharp and effective hook



R57040

## Pudendal Needles I R57040

Rocket® Pudendal Needles provide a safe and controllable method of delivering a Pudendal nerve block. The needle is disposable, has a luer-lock socket and is delivered sterile.

- Disposable
- Sterile
- Luer-lock socket



## Holder Scalpel I 9603

Made of stainless steel. Designed for easy handling and cleaning/autoclaving.

- Reusable
- Made of stainless steel
- Autoclavable

## Scalpel I 9602

Scalpel for precise cuts when performing blood sampling tests on the scalp of the baby. Made of stainless steel and environmentally friendly plastic.

- Sterile
- With cut-protection (remove before use)



## Holder Capillary Tube I 9605

Designed to hold the Capillary Tube 9604. Both ends can be used.

- Sterile

## Capillary Tube I 9604

Na-Heparinized tube for blood sampling test on baby's scalp. Used together with Holder for Capillary Tube 9605.

- Sterile
- 3 tubes/pouch
- Volume 75 ul



9607

9607S

## Cord Ring I 9607/9607S

Simple and easy-to-use umbilical cord ring which is available in Latex or Latex-free (9607S)

- Sterile
- Very elastic latex
- Available in Latex-free (9607S)



303100

## Cord Clamp I 303100

Cord Clamp for a safe and permanent clamping of the umbilical cord. Prevents loss of blood and shrivelled umbilical cord. Safe locking system. Keeps the cord firmly in the same position and does not slip of.

- Sterile
- Easy to use



434009

## Universal pump I 434009

Universal pump 434009 is advantageously used for large pack or dunk of ultrasonic gel, slider or similar gels and creams.

- Accessories for 5 and 1 kg cans
- It meets the requirements for substances that come into contact with food
- Very durable pump



434010

## Universal Dispenser I 434010

Universal dispenser 434010 is advantageously used for large pack or dunk of ultrasonic gel, slider or similar gels and creams.

- A light hand printing provides the desired amount
- It meets the requirements for materials that come in contact with food
- Hygienic and without spillage
- The contents of the packaging are utilized to the last drop and no cumbersome bottle refills are avoided



## Tube specula I 1340/1341

Wing Plast Tube Specula is a vaginal instrument in tube design, manufactured in Eco-friendly disposable plastic. The instrument are used to take blood tests on a child's scalp before delivery.

- Stable construction
- Smooth surfaces and smoothly rounded edges
- Smooth surface which facilitates insertion
- Easy to handle



## Aminoscope I 1342/1343/1344

Amnioscope is a flexible disposable instrument used in a variety of applications surveys and samples. The amnioscope series includes sizes; Large (long) and Large (short). Large (short) is the one with the highest consumption.

- Easy to handle
- Low friction surface for easy insertion
- Smooth surfaces and smoothly rounded edges
- Proven quality and strength
- Easy grip
- The battery is easy to separate for recycling



4171



4161



4151

## Ultrasound gel

4150/4151/4155/4160/4161/4170/4171

A specially formulated cosmetic based gel recommended for all procedures where a viscous gel is required. Available in sizes 260 g, 1 kg and 5 kg. Also available as bluecolored. And in Sterile packaging 20 g.

- Contains no salt
- Non-irritating - even after many scans
- Excellent viscosity (80,000 RPM5-STD 18C)
- Fast transfer rate
- Perfect acoustic impedance
- Extremely effective

Wing Plast offers a wide range of obstetric products, our easy to use and accurate Pregnancy calendar, Pregdate measurer and Cord rings are highly appreciated by midwives and healthcare professionals. Our Ultrasound gel is an extremely effective contact gel that is available in different sizes and has excellent viscosity (80,000 RPM5-STD 18C). Gynae Protection Sheet, Amniotic Membrane Perforator Cetrohook and Transducer bands are also available in our assortment. Take a look at our catalogue . If you have any questions about the products, please feel free to contact us.

**If you have an innovative idea in healthcare and want help with development or manufacturing, please contact us, we'll help you get your idea into reality.**



Wing Plast AB

**Wing Plast** - *Quality Medical Instruments*

Phone +46 (0)371 330 30, Fax +46 (0)371 330 31 , e-mail: [info@wingplast.se](mailto:info@wingplast.se), [www.wingplast.se](http://www.wingplast.se)