

ORIPLAST®



Surgical Suction
RELIABLE CLARIFICATION

Surgical Suction

FLEXIBLE SURGICAL SUCTION SYSTEMS

Systems for any need

Surgical suction systems are designed for the suction of fluids and tissue residues. They are mainly used during surgery to keep the operating area clean and maintain an overview or remove rinsing fluid that has been applied.


Surgical suction systems remove solid and liquid substances quickly and reliably without clogging.


Our connection tubes offer numerous connection options with a **wide variety of adapters**. The following connectors can be found in our range: Sherwood, funnel, silicone and fingertip. There is a choice between one-sided and two-sided connectors.

Adapter

The protection cap allows **uncomplicated and fast closure** of the respective adapter or connector. These also protect the user from contamination.

Protection cap





Suction tube


The length of the suction tubes varies from 1 to 5 m and allows optimal and reliable connection with a variety of adapters. In addition, these are available in different lumen sizes. The transparent disposable tubes are made of **high-quality, kink-resistant material** and guarantee a high flow rate

Fingertip

The fingertip can be used as a connection option. It is used for suction regulation and offers further connection options to the integrated step connector. Suction systems can be customised on request.

Suction Complete Set

The set can be used for an extraction of body fluids. It consists of a suction tube, a Redon drain and a disposable Redon needle.



Quad Pack

Tube: with slide clamp

Redon drain: 500 mm, 140 mm crossperforated, with an X-ray contrast line and length marking

Redon needle: straight form



Art.-No.	PCN	eClass	Tube length (mm)	Secretion stop	Connector	FG	PU
VK00512	17486277	34-24-13-90	1080	•	cut down connector	10	50

Suction Tube

Surgical suction tubes are suitable for the suction of body fluids. These are available in a different range of design and size.



Suction Connecting Tube


- for connecting suction pump and handpiece



Art.-No.	PCN	eClass	FG	Tube length (mm)	Connector	Endcap (without thread)	PU
VK00005	17483273	34-24-13-05	18	1800	two-sided silicone connector	-	100
VK00006	17483296	34-24-13-05	18	1800	silicone connector, LL pos	-	50
VK00129 ¹	17484255	34-24-13-05	25	1500	two-sided Sherwood	-	50
VK00128	17484249	34-24-13-05	25	1600	two-sided silicone connector	-	50
VK00131	17484278	34-24-13-05	25	1600	silicone connector, fingertip	-	50
VK00527	17486366	34-24-13-05	25	1600	one-side Sherwood	-	50
VK00529	17486372	34-24-13-05	25	1600	Sherwood, fingertip	-	50
VK00516	17486308	34-24-13-05	25	2000	two-sided Sherwood	-	50
VK00112	17484166	34-24-13-05	25	2100	two-sided funnel connector (grey)	-	50
VK00518	17486314	34-24-13-05	25	2400	two-sided Sherwood	-	50
VK00113	17484172	34-24-13-05	25	3000	two-sided funnel connector (grey)	-	50
VK00119	17484195	34-24-13-05	25	3000	funnel connector (grey)	-	50
VK00519	17486320	34-24-13-05	25	3000	two-sided Sherwood	-	50
VK00523 ^{2/3}	17486337	34-24-13-05	25	3000	two-sided Sherwood	-	50
VK00524 ³	17486343	34-24-13-05	25	3000	Sherwood, silicone connector	-	50
VK00120	17484203	34-24-13-05	25	3500	funnel connector (grey), Sherwood	-	50

Art.-No.	PCN	eClass	FG	Tube length (mm)	Connector	Endcap (without thread)	PU
VK00007	17483304	34-24-13-05	27	1000	two-sided silicone connector	-	50
VK00127	17484232	34-24-13-05	27	1300	two-sided silicone connector	-	50
VK00126	17484226	34-24-13-05	27	1900	two-sided silicone connector	-	50
VK00130	17484261	34-24-13-05	27	3000	silicone connector, fingertip	-	50
VK00001 ³	17483267	34-24-13-05	33	1800	silicone connector, fingertip	2x blue	50
VK00111	17484143	34-24-13-05	40	1900	step connector	-	25

1 with slide clamp 2 soft tube 3 DEHP-free




Connection Tube

Tube: soft material



Art.-No.	PCN	eClass	FG	Tube length (mm)	PU
VK01060	17489413	34-29-01-91	25	2200	50
VK00115	17484189	34-24-13-05	25	3000	50



Suction Tube

- with Y-connector



Art.-No.	PCN	eClass	FG	Tube length (mm)	PU
VK00916	17488537	34-29-01-01	28	280	50

Fingertip

The fingertip is an essential aid for secretion suction. The round tip part makes it possible to open or close the fingertip, which ultimately regulates and guarantees gentle suction.



Fingertip

- with round thumb-control valve




STERILE | EO

Art.-No.	PCN	eClass	Connector	PU
VK00950	17488796	34-24-13-07	-	100
VK00951	17488804	34-25-11-92	LL pos	100
VK00952	17488810	34-24-13-07	Luer-Steck	100

Secretion Bag

The latex-free collection bag is used to collect secretion and, with the help of a bottom drain valve, makes it easy and fast to take samples.



OP Secretion Bag

- tube length 900 mm
- volume 2000 ml
- with non-return valve
- bottom drain valve



STERILE | EO



Art.-No.	PCN	eClass	Connector	PU
VK00893	17488371	34-24-13-91	Luer-Steck with Polyisopren (Set), adapter	100

Suction

The Yankauer suction device provides a gentle and efficient suction method without aspiration of surrounding tissue. It can be connected to a fingertip, which allows a regulation of vacuum by thumb or finger pressure.




Yankauer

- angled suction tip
- without suction hose



Art.-No.	PCN	eClass	FG	Length (mm)	PU
VK00228	17484982	34-24-13-01	12	275	100
VK00227	17484976	34-24-13-01	18	280	100
VK00226	17484953	34-24-13-01	22	275	100

Adapter



Tube Adapter



Art.-No.	PCN	eClass	FG	Length (mm)	Connector	PU
VK00515	17486283	34-24-13-92	25	150	two-sided Sherwood	100