

CIMPAX® Fluid trap



Fluid trap

CIMPAX® Fluid Trap acts as a protective barrier, preventing the passage of fluids into the main filter of the smoke evacuation system and thereby maintaining the integrity of the filtration system.

CIMPAX® Fluid Trap reduces maintenance costs and extends the lifespan of the smoke evacuation equipment, contributing to a safer and more effective surgical smoke free environment.

- Translucent and lightweight design
- Small and compact - less waste material
- Clear measurement increments (holds 100ml)
- Easy to attach and remove
- Multiple connection options
- Integrated protection caps for safe disposal
- Single use
- Packaging: Box of 20 pcs. NON STERILE
- Ref. code 100ml: C32-30-010